

产品名称: FOXA2 Rabbit Monoclonal Antibody
产品货号: AMRe86635



产品概述 (Summary)

产品名称 (Production Name)	FOXA2 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC/IF
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	FOXA2
别名 (Alternative Names)	HNF3B; TCF3B
基因 ID (Gene ID)	3170
蛋白 ID (SwissProt ID)	Q9Y261.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000, ICC/IF 1:200-1:500
蛋白分子量 (Molecular Weight)	Calculated MW:48 kDa; Observed MW:50 kDa

研究背景 (Background)


产品名称: FOXA2 Rabbit Monoclonal Antibody
产品货号: AMRe86635



This gene encodes a member of the forkhead class of DNA-binding proteins. These hepatocyte nuclear factors are transcriptional activators for liver-specific genes such as albumin and transthyretin, and they also interact with chromatin. Similar family members in mice have roles in the regulation of metabolism and in the differentiation of the pancreas and liver. This gene has been linked to sporadic cases of maturity-onset diabetes of the young. Transcript variants encoding different isoforms have been identified for this gene. [provided by RefSeq, Oct 2008]

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of extracts from SW480 cells using AMRe86635 at 1:1000.

注意事项 (Note)

For research use only .