

产品名称: Baf180 Rabbit Monoclonal Antibody
产品货号: AMRe86640



产品概述 (Summary)

产品名称 (Production Name)	Baf180 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	IHC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Baf180
别名 (Alternative Names)	PB1; BAF180
基因 ID (Gene ID)	55193
蛋白 ID (SwissProt ID)	Q86U86.

产品应用 (Application)

稀释比 (Dilution Ratio)	IHC 1:500-1:2000
蛋白分子量 (Molecular Weight)	Calculated MW:193 kDa; Observed MW:193 kDa

研究背景 (Background)


产品名称: Baf180 Rabbit Monoclonal Antibody
产品货号: AMRe86640



This locus encodes a subunit of ATP-dependent chromatin-remodeling complexes. The encoded protein has been identified as an integral component of complexes necessary for ligand-dependent transcriptional activation by nuclear hormone receptors. Mutations at this locus have been associated with primary clear cell renal cell carcinoma. [provided by RefSeq, Feb 2012]

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of extracts from HeLa cells using AMRe86640 at 1:1000.

注意事项 (Note)

For research use only .