

产品名称: BCL-6 Rabbit Monoclonal Antibody
产品货号: AMRe86645



产品概述 (Summary)

产品名称 (Production Name)	BCL-6 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	BCL-6
别名 (Alternative Names)	BCL5; LAZ3; BCL6A; ZNF51; ZBTB27
基因 ID (Gene ID)	604
蛋白 ID (SwissProt ID)	P41182.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:5000,IHC 1:200-1:2000
蛋白分子量 (Molecular Weight)	Calculated MW:78 kDa; Observed MW:78 kDa

研究背景 (Background)


产品名称: BCL-6 Rabbit Monoclonal Antibody
产品货号: AMRe86645



The protein encoded by this gene is a zinc finger transcription factor and contains an N-terminal POZ domain. This protein acts as a sequence-specific repressor of transcription, and has been shown to modulate the transcription of STAT-dependent IL-4 responses of B cells. This protein can interact with a variety of POZ-containing proteins that function as transcription corepressors. This gene is found to be frequently translocated and hypermutated in diffuse large-cell lymphoma (DLCL), and may be involved in the pathogenesis of DLCL. Alternatively spliced transcript variants encoding different protein isoforms have been found for this gene. [provided by RefSeq, Aug 2015]

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of extracts from Raji cells using AMRe86645 at 1:1000.

注意事项 (Note)

For research use only .