

产品名称: HOXB13 Rabbit Monoclonal Antibody
产品货号: AMRe86646



产品概述 (Summary)

产品名称 (Production Name)	HOXB13 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,IP
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	HOXB13
别名 (Alternative Names)	PSGD
基因 ID (Gene ID)	10481
蛋白 ID (SwissProt ID)	Q92826.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:5000,IHC 1:500-1:2000,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW:31 kDa; Observed MW:31 kDa

研究背景 (Background)


产品名称: HOXB13 Rabbit Monoclonal Antibody
产品货号: AMRe86646



This gene encodes a transcription factor that belongs to the homeobox gene family. Genes of this family are highly conserved among vertebrates and essential for vertebrate embryonic development. This gene has been implicated to play a role in fetal skin development and cutaneous regeneration. In mice, a similar gene was shown to exhibit temporal and spatial colinearity in the main body axis of the embryo, but was not expressed in the secondary axes, which suggests functions in body patterning along the axis. This gene and other HOXB genes form a gene cluster at chromosome the 17q21-22 region. [provided by RefSeq, Jul 2008]

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of extracts from LNCaP cells using AMRe86646 at 1:1000.

注意事项 (Note)

For research use only .