

产品名称: Somatostatin Receptor 2 Rabbit Monoclonal Antibody
产品货号: AMRe86654



产品概述 (Summary)

产品名称 (Production Name)	Somatostatin Receptor 2 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Somatostatin Receptor 2
别名 (Alternative Names)	SSTR2
基因 ID (Gene ID)	6752
蛋白 ID (SwissProt ID)	P30874.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:5000,IHC 1:100-1:1000
蛋白分子量 (Molecular Weight)	Calculated MW:41 kDa; Observed MW:60 kDa

产品名称: Somatostatin Receptor 2 Rabbit Monoclonal Antibody
产品货号: AMRe86654




研究背景 (Background)

Somatostatin acts at many sites to inhibit the release of many hormones and other secretory proteins. The biologic effects of somatostatin are probably mediated by a family of G protein-coupled receptors that are expressed in a tissue-specific manner. SSTR2 is a member of the superfamily of receptors having seven transmembrane segments and is expressed in highest levels in cerebrum and kidney. [provided by RefSeq, Jul 2008]

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of extracts from Mouse brain tissue using AMRe86654 at 1:1000.

注意事项 (Note)

For research use only .