

产品名称: HMGCS2 Rabbit Monoclonal Antibody
产品货号: AMRe86662



产品概述 (Summary)

产品名称 (Production Name)	HMGCS2 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	HMGCS2
别名 (Alternative Names)	3-hydroxy-3-methylglutaryl coenzyme A synthase
基因 ID (Gene ID)	3158
蛋白 ID (SwissProt ID)	P54868.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:5000,IHC 1:50-1:200
蛋白分子量 (Molecular Weight)	Calculated MW:57 kDa; Observed MW:57 kDa

研究背景 (Background)


产品名称: HMGCS2 Rabbit Monoclonal Antibody
产品货号: AMRe86662



The protein encoded by this gene belongs to the HMG-CoA synthase family. It is a mitochondrial enzyme that catalyzes the first reaction of ketogenesis, a metabolic pathway that provides lipid-derived energy for various organs during times of carbohydrate deprivation, such as fasting. Mutations in this gene are associated with HMG-CoA synthase deficiency. Alternatively spliced transcript variants encoding different isoforms have been found for this gene.[provided by RefSeq, Oct 2009]

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of extracts from Caco-2 cells using AMRe86662 at 1:1000.

注意事项 (Note)

For research use only .