

产品名称: PSMA Rabbit Monoclonal Antibody
产品货号: AMRe86671



产品概述 (Summary)

产品名称 (Production Name)	PSMA Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,IP
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	PSMA
别名 (Alternative Names)	PSM; FGCP; FOLH; GCP2; PSMA; mGCP; GCPII; NAALAD1; NAALAdase
基因 ID (Gene ID)	2346
蛋白 ID (SwissProt ID)	Q04609.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:5000-1:50000,IHC 1:200-1:500,IP 1:10-1:100
蛋白分子量 (Molecular Weight)	Calculated MW:84 kDa; Observed MW:100 kDa

研究背景 (Background)


产品名称: PSMA Rabbit Monoclonal Antibody
产品货号: AMRe86671



This gene encodes a type II transmembrane glycoprotein belonging to the M28 peptidase family. The protein acts as a glutamate carboxypeptidase on different alternative substrates, including the nutrient folate and the neuropeptide N-acetyl-L-aspartyl-L-glutamate and is expressed in a number of tissues such as prostate, central and peripheral nervous system and kidney. A mutation in this gene may be associated with impaired intestinal absorption of dietary folates, resulting in low blood folate levels and consequent hyperhomocysteinemia. Expression of this protein in the brain may be involved in a number of pathological conditions associated with glutamate excitotoxicity. In the prostate the protein is up-regulated in cancerous cells and is used as an effective diagnostic and prognostic indicator of prostate cancer. This gene likely arose from a duplication event of a nearby chromosomal region. Alternative splicing gives rise to multiple transcript variants encoding several different isoforms. [provided by RefSeq, Jul 2010]

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of extracts from LNCaP cells using AMRe86671 at 1:1000.

注意事项 (Note)

For research use only .