

产品名称: FGR Rabbit Monoclonal Antibody  
产品货号: AMRe86681



## 产品概述 (Summary)

|                        |                                        |
|------------------------|----------------------------------------|
| 产品名称 (Production Name) | FGR Rabbit Monoclonal Antibody         |
| 描述 (Description)       | Recombinant rabbit monoclonal antibody |
| 宿主 (Host)              | Rabbit                                 |
| 应用 (Application)       | WB                                     |
| 种属反应性 (Reactivity)     | Human                                  |

## 产品性能 (Performance)

|                     |                                                                                                                                                              |
|---------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 偶联物 (Conjugation)   | Unconjugated                                                                                                                                                 |
| 修饰 (Modification)   | Unmodified                                                                                                                                                   |
| 同种型 (Isotype)       | IgG                                                                                                                                                          |
| 克隆 (Clonality)      | Monoclonal                                                                                                                                                   |
| 形式 (Form)           | Liquid                                                                                                                                                       |
| 存放说明 (Storage)      | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.                                                                     |
| 储存溶液 (Buffer)       | Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt. |
| 纯化方式 (Purification) | Affinity Purification                                                                                                                                        |

## 免疫原信息 (Immunogen)

|                        |                                                           |
|------------------------|-----------------------------------------------------------|
| 基因名 (Gene Name)        | FGR                                                       |
| 别名 (Alternative Names) | SRC2; c-fgr; c-src2; p55-Fgr; p58-Fgr; p55c-fgr; p58c-fgr |
| 基因 ID (Gene ID)        | 2268                                                      |
| 蛋白 ID (SwissProt ID)   | P09769.                                                   |

## 产品应用 (Application)

|                          |                                          |
|--------------------------|------------------------------------------|
| 稀释比 (Dilution Ratio)     | WB 1:1000-1:5000                         |
| 蛋白分子量 (Molecular Weight) | Calculated MW:60 kDa; Observed MW:60 kDa |

## 研究背景 (Background)


**产品名称: FGR Rabbit Monoclonal Antibody**  
**产品货号: AMRe86681**



This gene is a member of the Src family of protein tyrosine kinases (PTKs). The encoded protein contains N-terminal sites for myristylation and palmitoylation, a PTK domain, and SH2 and SH3 domains which are involved in mediating protein-protein interactions with phosphotyrosine-containing and proline-rich motifs, respectively. The protein localizes to plasma membrane ruffles, and functions as a negative regulator of cell migration and adhesion triggered by the beta-2 integrin signal transduction pathway. Infection with Epstein-Barr virus results in the overexpression of this gene. Multiple alternatively spliced variants, encoding the same protein, have been identified. [provided by RefSeq, Jul 2008]

## 研究领域 (Research Area)

## 图片 (Image Data)



Western blot analysis of extracts from Raji cells using AMRe86681 at 1:1000.

## 注意事项 (Note)

For research use only .