

**产品名称: Glutathione S Transferase kappa 1 Rabbit
Monoclonal Antibody**
产品货号: AMRe86695



产品概述 (Summary)

产品名称 (Production Name)	Glutathione S Transferase kappa 1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,FC,IP
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Glutathione S Transferase kappa 1
别名 (Alternative Names)	GST; GST13; hGSTK1; GSTK1-1; GST13-13; GST 13-13
基因 ID (Gene ID)	373156
蛋白 ID (SwissProt ID)	Q9Y2Q3.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:100-1:200,ICC/IF 1:20-1:50,FC 1:20-1:50,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW:26 kDa; Observed MW:26 kDa

**产品名称: Glutathione S Transferase kappa 1 Rabbit
Monoclonal Antibody**
产品货号: AMRe86695




研究背景 (Background)

This gene encodes a member of the kappa class of the glutathione transferase superfamily of enzymes that function in cellular detoxification. The encoded protein is localized to the peroxisome and catalyzes the conjugation of glutathione to a wide range of hydrophobic substrates facilitating the removal of these compounds from cells. Alternative splicing results in multiple transcript variants.[provided by RefSeq, Jan 2009]

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of extracts from Raji cells using AMRe86695 at 1:1000.

注意事项 (Note)

For research use only .