

产品名称: PON3 Rabbit Monoclonal Antibody  
产品货号: AMRe86706



## 产品概述 (Summary)

|                        |                                        |
|------------------------|----------------------------------------|
| 产品名称 (Production Name) | PON3 Rabbit Monoclonal Antibody        |
| 描述 (Description)       | Recombinant rabbit monoclonal antibody |
| 宿主 (Host)              | Rabbit                                 |
| 应用 (Application)       | WB                                     |
| 种属反应性 (Reactivity)     | Human, Mouse                           |

## 产品性能 (Performance)

|                     |                                                                                                                                                              |
|---------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 偶联物 (Conjugation)   | Unconjugated                                                                                                                                                 |
| 修饰 (Modification)   | Unmodified                                                                                                                                                   |
| 同种型 (Isotype)       | IgG                                                                                                                                                          |
| 克隆 (Clonality)      | Monoclonal                                                                                                                                                   |
| 形式 (Form)           | Liquid                                                                                                                                                       |
| 存放说明 (Storage)      | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.                                                                     |
| 储存溶液 (Buffer)       | Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt. |
| 纯化方式 (Purification) | Affinity Purification                                                                                                                                        |

## 免疫原信息 (Immunogen)

|                        |                               |
|------------------------|-------------------------------|
| 基因名 (Gene Name)        | PON3                          |
| 别名 (Alternative Names) | Serum paraoxonase/lactonase 3 |
| 基因 ID (Gene ID)        | 5446                          |
| 蛋白 ID (SwissProt ID)   | Q15166.                       |

## 产品应用 (Application)

|                          |                                          |
|--------------------------|------------------------------------------|
| 稀释比 (Dilution Ratio)     | WB 1:1000-1:5000                         |
| 蛋白分子量 (Molecular Weight) | Calculated MW:40 kDa; Observed MW:40 kDa |

## 研究背景 (Background)


**产品名称: PON3 Rabbit Monoclonal Antibody**  
**产品货号: AMRe86706**



This gene is a member of the paraoxonase family and lies in a cluster on chromosome 7 with the other two family members. The encoded protein is secreted into the bloodstream and associates with high-density lipoprotein (HDL). The protein also rapidly hydrolyzes lactones and can inhibit the oxidation of low-density lipoprotein (LDL), a function that is believed to slow the initiation and progression of atherosclerosis. Alternatively spliced variants which encode different protein isoforms have been described; however, only one has been fully characterized. [provided by RefSeq, Jul 2008]

## 研究领域 (Research Area)

## 图片 (Image Data)



Western blot analysis of extracts from HepG2 cells using AMRe86706 at 1:1000.

## 注意事项 (Note)

For research use only .