

产品名称: LC3B Rabbit Monoclonal Antibody
产品货号: AMRe86710



产品概述 (Summary)

产品名称 (Production Name)	LC3B Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF
种属反应性 (Reactivity)	Human, Mouse, Rat, Other

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	LC3B
别名 (Alternative Names)	LC3B; ATG8F; MAP1LC3B-a; MAP1A/1BLC3
基因 ID (Gene ID)	81631, 67443, 64862
蛋白 ID (SwissProt ID)	Q9GZQ8, Q9CQV6, Q62625.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:5000,IHC 1:500-1:2000,ICC/IF 1:500-1:2000
蛋白分子量 (Molecular Weight)	Calculated MW:15 kDa; Observed MW:14,16 kDa

产品名称: LC3B Rabbit Monoclonal Antibody
产品货号: AMRe86710




研究背景 (Background)


The product of this gene is a subunit of neuronal microtubule-associated MAP1A and MAP1B proteins, which are involved in microtubule assembly and important for neurogenesis. Studies on the rat homolog implicate a role for this gene in autophagy, a process that involves the bulk degradation of cytoplasmic component. [provided by RefSeq, Jul 2008]

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of extracts from U-87 MG cells using AMRe86710 at 1:1000.



Western blot analysis of extracts from MCF-7 cells and Rat muscle tissue using AMRe86710 at 1:4000.

注意事项 (Note)

For research use only .