

产品名称: Calcineurin A Rabbit Monoclonal Antibody
产品货号: AMRe86712



产品概述 (Summary)

产品名称 (Production Name)	Calcineurin A Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC/IF, FC
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Calcineurin A
别名 (Alternative Names)	CALN; CCN1; CNA1; CALNA; PPP2B; CALNA1
基因 ID (Gene ID)	5530
蛋白 ID (SwissProt ID)	Q08209.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:2000-1:20000, ICC/IF 1:100-1:200, FC 1:10-1:100
蛋白分子量 (Molecular Weight)	Calculated MW:59 kDa; Observed MW:59 kDa

研究背景 (Background)


产品名称: Calcineurin A Rabbit Monoclonal Antibody
产品货号: AMRe86712



Calcium-dependent, calmodulin-stimulated protein phosphatase. Many of the substrates contain a PxIxIT motif. This subunit may have a role in the calmodulin activation of calcineurin. Dephosphorylates DNM1L, HSPB1 and SSH1.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of extracts from Mouse brain tissue using AMRe86712 at 1:1000.

注意事项 (Note)

For research use only .