

产品名称: DKK1 Rabbit Monoclonal Antibody
产品货号: AMRe86718



产品概述 (Summary)

产品名称 (Production Name)	DKK1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	DKK1
别名 (Alternative Names)	SK; DKK-1
基因 ID (Gene ID)	22943
蛋白 ID (SwissProt ID)	O94907.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:5000,IHC 1:200-1:1000
蛋白分子量 (Molecular Weight)	Calculated MW:29 kDa; Observed MW:28-40 kDa

研究背景 (Background)


产品名称: DKK1 Rabbit Monoclonal Antibody
产品货号: AMRe86718



This gene encodes a member of the dickkopf family of proteins. Members of this family are secreted proteins characterized by two cysteine-rich domains that mediate protein-protein interactions. The encoded protein binds to the LRP6 co-receptor and inhibits beta-catenin-dependent Wnt signaling. This gene plays a role in embryonic development and may be important in bone formation in adults. Elevated expression of this gene has been observed in numerous human cancers and this protein may promote proliferation, invasion and growth in cancer cell lines. [provided by RefSeq, Sep 2017]

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of extracts from PC-3 cells using AMRe86718 at 1:1000.

注意事项 (Note)

For research use only .