

产品名称: TEAD-1 Rabbit Monoclonal Antibody
产品货号: AMRe86731



产品概述 (Summary)

产品名称 (Production Name)	TEAD-1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,IP
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	TEAD-1
别名 (Alternative Names)	AA; REF1; TCF13; TEF-1; NTEF-1; TCF-13; TEAD-1
基因 ID (Gene ID)	7003
蛋白 ID (SwissProt ID)	P28347.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:5000,IHC 1:200-1:1000,IP 1:10-1:100
蛋白分子量 (Molecular Weight)	Calculated MW:48 kDa; Observed MW:48 kDa

研究背景 (Background)


产品名称: TEAD-1 Rabbit Monoclonal Antibody
产品货号: AMRe86731



This gene encodes a ubiquitous transcriptional enhancer factor that is a member of the TEA/ATTS domain family. This protein directs the transactivation of a wide variety of genes and, in placental cells, also acts as a transcriptional repressor. Mutations in this gene cause Sveinsson's chorioretinal atrophy. Additional transcript variants have been described but their full-length natures have not been experimentally verified. [provided by RefSeq, May 2010]

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of extracts from A549, HEK293T, HeLa, H9C2 cells using AMRe86731 at 1:1000.

注意事项 (Note)

For research use only .