

**产品名称: HRPT2 Rabbit Monoclonal Antibody**  
**产品货号: AMRe86737**



## 产品概述 (Summary)

|                        |                                        |
|------------------------|----------------------------------------|
| 产品名称 (Production Name) | HRPT2 Rabbit Monoclonal Antibody       |
| 描述 (Description)       | Recombinant rabbit monoclonal antibody |
| 宿主 (Host)              | Rabbit                                 |
| 应用 (Application)       | WB,IHC,FC                              |
| 种属反应性 (Reactivity)     | Human,Mouse,Rat                        |

## 产品性能 (Performance)

|                     |                                                                                                                                                              |
|---------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 偶联物 (Conjugation)   | Unconjugated                                                                                                                                                 |
| 修饰 (Modification)   | Unmodified                                                                                                                                                   |
| 同种型 (Isotype)       | IgG                                                                                                                                                          |
| 克隆 (Clonality)      | Monoclonal                                                                                                                                                   |
| 形式 (Form)           | Liquid                                                                                                                                                       |
| 存放说明 (Storage)      | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.                                                                     |
| 储存溶液 (Buffer)       | Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt. |
| 纯化方式 (Purification) | Affinity Purification                                                                                                                                        |

## 免疫原信息 (Immunogen)

|                        |                                         |
|------------------------|-----------------------------------------|
| 基因名 (Gene Name)        | HRPT2                                   |
| 别名 (Alternative Names) | HYX; FIHP; HPTJT; HRPT1; HRPT2; C1orf28 |
| 基因 ID (Gene ID)        | 79577                                   |
| 蛋白 ID (SwissProt ID)   | Q6P1J9.                                 |

## 产品应用 (Application)

|                          |                                                 |
|--------------------------|-------------------------------------------------|
| 稀释比 (Dilution Ratio)     | WB 1:500-1:2000,IHC 1:200-1:1000,FC 1:200-1:500 |
| 蛋白分子量 (Molecular Weight) | Calculated MW:61 kDa; Observed MW:61 kDa        |

## 研究背景 (Background)


**产品名称: HRPT2 Rabbit Monoclonal Antibody**  
**产品货号: AMRe86737**



This gene encodes a tumor suppressor that is involved in transcriptional and post-transcriptional control pathways. The protein is a component of the the PAF protein complex, which associates with the RNA polymerase II subunit POLR2A and with a histone methyltransferase complex. This protein appears to facilitate the association of 3' mRNA processing factors with actively-transcribed chromatin. Mutations in this gene have been linked to hyperparathyroidism-jaw tumor syndrome, familial isolated hyperparathyroidism, and parathyroid carcinoma. [provided by RefSeq, Jul 2009]

## 研究领域 (Research Area)

## 图片 (Image Data)



Western blot analysis of extracts from HepG2 cells using AMRe86737 at 1:1000.

## 注意事项 (Note)

For research use only .