

产品名称: p73 Rabbit Monoclonal Antibody  
产品货号: AMRe86739



## 产品概述 (Summary)

|                        |                                        |
|------------------------|----------------------------------------|
| 产品名称 (Production Name) | p73 Rabbit Monoclonal Antibody         |
| 描述 (Description)       | Recombinant rabbit monoclonal antibody |
| 宿主 (Host)              | Rabbit                                 |
| 应用 (Application)       | WB,IHC,IP                              |
| 种属反应性 (Reactivity)     | Human                                  |

## 产品性能 (Performance)

|                     |                                                                                                                                                              |
|---------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 偶联物 (Conjugation)   | Unconjugated                                                                                                                                                 |
| 修饰 (Modification)   | Unmodified                                                                                                                                                   |
| 同种型 (Isotype)       | IgG                                                                                                                                                          |
| 克隆 (Clonality)      | Monoclonal                                                                                                                                                   |
| 形式 (Form)           | Liquid                                                                                                                                                       |
| 存放说明 (Storage)      | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.                                                                     |
| 储存溶液 (Buffer)       | Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt. |
| 纯化方式 (Purification) | Affinity Purification                                                                                                                                        |

## 免疫原信息 (Immunogen)

|                        |         |
|------------------------|---------|
| 基因名 (Gene Name)        | p73     |
| 别名 (Alternative Names) | P73     |
| 基因 ID (Gene ID)        | 7161    |
| 蛋白 ID (SwissProt ID)   | O15350. |

## 产品应用 (Application)

|                          |                                              |
|--------------------------|----------------------------------------------|
| 稀释比 (Dilution Ratio)     | WB 1:500-1:2000,IHC 1:200-1:500,IP 1:20-1:50 |
| 蛋白分子量 (Molecular Weight) | Calculated MW:70 kDa; Observed MW:80,70 kDa  |

## 研究背景 (Background)


**产品名称: p73 Rabbit Monoclonal Antibody**  
**产品货号: AMRe86739**



This gene encodes a member of the p53 family of transcription factors involved in cellular responses to stress and development. It maps to a region on chromosome 1p36 that is frequently deleted in neuroblastoma and other tumors, and thought to contain multiple tumor suppressor genes. The demonstration that this gene is monoallelically expressed (likely from the maternal allele), supports the notion that it is a candidate gene for neuroblastoma. Many transcript variants resulting from alternative splicing and/or use of alternate promoters have been found for this gene, but the biological validity and the full-length nature of some variants have not been determined. [provided by RefSeq, Feb 2011]

## 研究领域 (Research Area)

## 图片 (Image Data)



Western blot analysis of extracts from HeLa cells using AMRe86739 at 1:1000.

## 注意事项 (Note)

For research use only .