

产品名称: Smoothelin Rabbit Monoclonal Antibody
产品货号: AMRe86750



产品概述 (Summary)

产品名称 (Production Name)	Smoothelin Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Smoothelin
别名 (Alternative Names)	SMTN; SMSMO
基因 ID (Gene ID)	6525
蛋白 ID (SwissProt ID)	P53814.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:1000-1:5000
蛋白分子量 (Molecular Weight)	Calculated MW:99 kDa; Observed MW:110,59 kDa

研究背景 (Background)


产品名称: Smoothelin Rabbit Monoclonal Antibody
产品货号: AMRe86750



This gene encodes a structural protein that is found exclusively in contractile smooth muscle cells. It associates with stress fibers and constitutes part of the cytoskeleton. This gene is localized to chromosome 22q12.3, distal to the TUPLE1 locus and outside the DiGeorge syndrome deletion. Alternative splicing of this gene results in multiple transcript variants encoding distinct isoforms. [provided by RefSeq, May 2011]

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of extracts from Mouse colon tissue using AMRe86750 at 1:1000.

注意事项 (Note)

For research use only .