

产品名称: LKB1 Rabbit Monoclonal Antibody
产品货号: AMRe86759



产品概述 (Summary)

产品名称 (Production Name)	LKB1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IP
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	LKB1
别名 (Alternative Names)	PJS; LKB1; hLKB1
基因 ID (Gene ID)	6794
蛋白 ID (SwissProt ID)	Q15831.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IP 1:50-1:200
蛋白分子量 (Molecular Weight)	Calculated MW:49 kDa; Observed MW:49 kDa

研究背景 (Background)


产品名称: LKB1 Rabbit Monoclonal Antibody
产品货号: AMRe86759



This gene, which encodes a member of the serine/threonine kinase family, regulates cell polarity and functions as a tumor suppressor. Mutations in this gene have been associated with Peutz-Jeghers syndrome, an autosomal dominant disorder characterized by the growth of polyps in the gastrointestinal tract, pigmented macules on the skin and mouth, and other neoplasms. Alternate transcriptional splice variants of this gene have been observed but have not been thoroughly characterized. [provided by RefSeq, Jul 2008]

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of extracts from K562 cells using AMRe86759 at 1:100.

注意事项 (Note)

For research use only .