

产品名称: CD4 Rabbit Monoclonal Antibody
产品货号: AMRe86776



产品概述 (Summary)

产品名称 (Production Name)	CD4 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,IP
种属反应性 (Reactivity)	Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	CD4
别名 (Alternative Names)	L3T4; Ly-4
基因 ID (Gene ID)	12504
蛋白 ID (SwissProt ID)	P06332.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:200-1:1000,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW:51 kDa; Observed MW:51 kDa

研究背景 (Background)


产品名称: CD4 Rabbit Monoclonal Antibody
产品货号: AMRe86776



Integral membrane glycoprotein that plays an essential role in the immune response and serves multiple functions in responses against both external and internal offenses. In T-cells, functions primarily as a coreceptor for MHC class II molecule:peptide complex. The antigens presented by class II peptides are derived from extracellular proteins while class I peptides are derived from cytosolic proteins. Interacts simultaneously with the T-cell receptor (TCR) and the MHC class II presented by antigen presenting cells (APCs). In turn, recruits the Src kinase LCK to the vicinity of the TCR-CD3 complex. LCK then initiates different intracellular signaling pathways by phosphorylating various substrates ultimately leading to lymphokine production, motility, adhesion and activation of T-helper cells.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of extracts from Mouse thymus tissue using AMRe86776 at 1:1000.

注意事项 (Note)

For research use only .