

产品名称: IGFBP4 Rabbit Monoclonal Antibody
产品货号: AMRe86780



产品概述 (Summary)

产品名称 (Production Name)	IGFBP4 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,FC,IP
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	IGFBP4
别名 (Alternative Names)	BP-4; IBP4; IGFBP-4; HT29-IGFBP
基因 ID (Gene ID)	3487
蛋白 ID (SwissProt ID)	P22692.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,FC 1:200-1:500,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW:28 kDa; Observed MW:28 kDa

研究背景 (Background)


产品名称: IGFBP4 Rabbit Monoclonal Antibody
产品货号: AMRe86780




This gene is a member of the insulin-like growth factor binding protein (IGFBP) family and encodes a protein with an IGFBP domain and a thyroglobulin type-I domain. The protein binds both insulin-like growth factors (IGFs) I and II and circulates in the plasma in both glycosylated and non-glycosylated forms. Binding of this protein prolongs the half-life of the IGFs and alters their interaction with cell surface receptors. [provided by RefSeq, Jul 2008]

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of extracts from HT-29 cells using AMRe86780 at 1:1000.



Western blot analysis of extracts from Rat heart, Rat muscle tissue using AMRe86780 at 1:1000.

注意事项 (Note)

For research use only .