

产品名称: ATF3 Rabbit Monoclonal Antibody
产品货号: AMRe86781



产品概述 (Summary)

产品名称 (Production Name)	ATF3 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC/IF, IP
种属反应性 (Reactivity)	Human, Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	ATF3
别名 (Alternative Names)	Activating transcription factor 3
基因 ID (Gene ID)	467
蛋白 ID (SwissProt ID)	P18847.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000, ICC/IF 1:100-1:200, IP 1:50-1:100
蛋白分子量 (Molecular Weight)	Calculated MW:21 kDa; Observed MW:21 kDa

研究背景 (Background)


产品名称: ATF3 Rabbit Monoclonal Antibody
产品货号: AMRe86781




This gene encodes a member of the mammalian activation transcription factor/cAMP responsive element-binding (CREB) protein family of transcription factors. This gene is induced by a variety of signals, including many of those encountered by cancer cells, and is involved in the complex process of cellular stress response. Multiple transcript variants encoding different isoforms have been found for this gene. It is possible that alternative splicing of this gene may be physiologically important in the regulation of target genes. [provided by RefSeq, Apr 2011]

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of extracts from HeLa, 293T cells using AMRe86781 at 1:1000.



Western blot analysis of extracts from RAW264.7 cells using AMRe86781 at 1:1000.

注意事项 (Note)

For research use only .