

产品名称: AMCase Rabbit Monoclonal Antibody
产品货号: AMRe86784



产品概述 (Summary)

产品名称 (Production Name)	AMCase Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	AMCase
别名 (Alternative Names)	CHIT2; AMCASE; TSA1902
基因 ID (Gene ID)	27159
蛋白 ID (SwissProt ID)	Q9BZP6.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:1000-1:5000
蛋白分子量 (Molecular Weight)	Calculated MW:52 kDa; Observed MW:52 kDa

研究背景 (Background)

产品名称: AMCase Rabbit Monoclonal Antibody
产品货号: AMRe86784




The protein encoded by this gene degrades chitin, which is found in the cell wall of most fungi as well as in arthropods and some nematodes. The encoded protein can also stimulate interleukin 13 expression, and variations in this gene can lead to asthma susceptibility. Several transcript variants encoding a few different isoforms have been found for this gene.

[provided by RefSeq, Apr 2012]

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of extracts from Mouse stomach tissue using AMRe86784 at 1:1000.

注意事项 (Note)

For research use only .