

产品名称: EGF Rabbit Monoclonal Antibody
产品货号: AMRe86789



产品概述 (Summary)

产品名称 (Production Name)	EGF Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	EGF
别名 (Alternative Names)	AI790464
基因 ID (Gene ID)	13645
蛋白 ID (SwissProt ID)	P01132.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:200-1:1000
蛋白分子量 (Molecular Weight)	Calculated MW:133 kDa; Observed MW:160 kDa

研究背景 (Background)


产品名称: EGF Rabbit Monoclonal Antibody
产品货号: AMRe86789



This gene encodes epidermal growth factor (EGF), the founding member of the EGF family of growth factors that are implicated in cell proliferation and differentiation. The encoded protein can localize to the membrane and function in juxtacrine signaling or undergo proteolytic processing to generate a soluble form of the hormone. Mice lacking the encoded protein do not exhibit an abnormal phenotype but transgenic mice overexpressing the encoded protein exhibit hypospermatogenesis. [provided by RefSeq, Jul 2016]

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of extracts from Mouse kidney tissue using AMRe86789 at 1:1000.

注意事项 (Note)

For research use only .