

产品名称: TLS/FUS Rabbit Monoclonal Antibody
产品货号: AMRe86791



产品概述 (Summary)

产品名称 (Production Name)	TLS/FUS Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,FC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	TLS/FUS
别名 (Alternative Names)	TLS; ALS6; ETM4; FUS1; POMP75; HNRNPP2
基因 ID (Gene ID)	2521
蛋白 ID (SwissProt ID)	P35637.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:200-1:1000,ICC/IF 1:100-1:200,FC 1:50-1:200
蛋白分子量 (Molecular Weight)	Calculated MW:53 kDa; Observed MW:70 kDa

研究背景 (Background)


产品名称: TLS/FUS Rabbit Monoclonal Antibody
产品货号: AMRe86791



This gene encodes a multifunctional protein component of the heterogeneous nuclear ribonucleoprotein (hnRNP) complex. The hnRNP complex is involved in pre-mRNA splicing and the export of fully processed mRNA to the cytoplasm. This protein belongs to the FET family of RNA-binding proteins which have been implicated in cellular processes that include regulation of gene expression, maintenance of genomic integrity and mRNA/microRNA processing. Alternative splicing results in multiple transcript variants. Defects in this gene result in amyotrophic lateral sclerosis type 6. [provided by RefSeq, Sep 2009]

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of extracts from HepG2 cells using AMRe86791 at 1:1000.

注意事项 (Note)

For research use only .