

产品名称: ODF1 Rabbit Monoclonal Antibody
产品货号: AMRe86810



产品概述 (Summary)

产品名称 (Production Name)	ODF1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	ODF1
别名 (Alternative Names)	Odfp
基因 ID (Gene ID)	18285
蛋白 ID (SwissProt ID)	Q61999.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:2000-1:20000
蛋白分子量 (Molecular Weight)	Calculated MW:28 kDa; Observed MW:28 kDa

研究背景 (Background)


产品名称: ODF1 Rabbit Monoclonal Antibody
产品货号: AMRe86810



Component of the outer dense fibers (ODF) of spermatozoa. ODF are filamentous structures located on the outside of the axoneme in the midpiece and principal piece of the mammalian sperm tail and may help to maintain the passive elastic structures and elastic recoil of the sperm tail.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of extracts from Mouse testis tissue using AMRe86810 at 1:1000.

注意事项 (Note)

For research use only .