

产品名称: RHOD Rabbit Monoclonal Antibody
产品货号: AMRe86811



产品概述 (Summary)

产品名称 (Production Name)	RHOD Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	RHOD
别名 (Alternative Names)	Rho; ARHD; RHOM; RHOHP1
基因 ID (Gene ID)	29984
蛋白 ID (SwissProt ID)	O00212.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000
蛋白分子量 (Molecular Weight)	Calculated MW:23 kDa; Observed MW:23 kDa

研究背景 (Background)


产品名称: RHOD Rabbit Monoclonal Antibody
产品货号: AMRe86811



Ras homolog, or Rho, proteins interact with protein kinases and may serve as targets for activated GTPase. They play a critical role in muscle differentiation. The protein encoded by this gene binds GTP and is a member of the small GTPase superfamily. It is involved in endosome dynamics and reorganization of the actin cytoskeleton, and it may coordinate membrane transport with the function of the cytoskeleton. Two transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jul 2014]

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of extracts from Mouse kidney tissue using AMRe86811 at 1:1000.

注意事项 (Note)

For research use only .