

产品名称: Nephrin Rabbit Monoclonal Antibody
产品货号: AMRe86830



产品概述 (Summary)

产品名称 (Production Name)	Nephrin Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,IP
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Nephrin
别名 (Alternative Names)	CNF; NPHN; nephrin
基因 ID (Gene ID)	4868
蛋白 ID (SwissProt ID)	O60500.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:200-1:1000,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW:135 kDa; Observed MW:200 kDa

研究背景 (Background)


产品名称: Nephrin Rabbit Monoclonal Antibody
产品货号: AMRe86830



This gene encodes a member of the immunoglobulin family of cell adhesion molecules that functions in the glomerular filtration barrier in the kidney. The gene is primarily expressed in renal tissues, and the protein is a type-1 transmembrane protein found at the slit diaphragm of glomerular podocytes. The slit diaphragm is thought to function as an ultrafilter to exclude albumin and other plasma macromolecules in the formation of urine. Mutations in this gene result in Finnish-type congenital nephrosis 1, characterized by severe proteinuria and loss of the slit diaphragm and foot processes.[provided by RefSeq, Oct 2009]

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of extracts from Mouse kidney tissue using AMRe86830 at 1:1000.

注意事项 (Note)

For research use only .