

产品名称: AGXT Rabbit Monoclonal Antibody
产品货号: AMRe86836



产品概述 (Summary)

产品名称 (Production Name)	AGXT Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,IP
种属反应性 (Reactivity)	Human,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	AGXT
别名 (Alternative Names)	AGT; PH1; SPT; AGT1; SPAT; TLH6; AGXT1
基因 ID (Gene ID)	189
蛋白 ID (SwissProt ID)	P21549.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:5000,IHC 1:100-1:200,IP 1:10-1:100
蛋白分子量 (Molecular Weight)	Calculated MW:43 kDa; Observed MW:43 kDa

研究背景 (Background)


产品名称: AGXT Rabbit Monoclonal Antibody
产品货号: AMRe86836



This gene is expressed only in the liver and the encoded protein is localized mostly in the peroxisomes, where it is involved in glyoxylate detoxification. Mutations in this gene, some of which alter subcellular targeting, have been associated with type I primary hyperoxaluria. [provided by RefSeq, Jul 2008]

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of extracts from HepG2 cells using AMRe86836 at 1:1000.

注意事项 (Note)

For research use only .