

产品名称: CD133 Rabbit Monoclonal Antibody
产品货号: AMRe86843



产品概述 (Summary)

产品名称 (Production Name)	CD133 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,FC,IP
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	CD133
别名 (Alternative Names)	RP41; AC133; CD133; MCDR2; STGD4; CORD12; PROML1; MSTP061
基因 ID (Gene ID)	8842
蛋白 ID (SwissProt ID)	O43490.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:200-1:1000,FC 1:50-1:200,IP 1:50-1:100
蛋白分子量 (Molecular Weight)	Calculated MW:97 kDa; Observed MW:133 kDa

研究背景 (Background)


产品名称: CD133 Rabbit Monoclonal Antibody
产品货号: AMRe86843



This gene encodes a pentaspan transmembrane glycoprotein. The protein localizes to membrane protrusions and is often expressed on adult stem cells, where it is thought to function in maintaining stem cell properties by suppressing differentiation. Mutations in this gene have been shown to result in retinitis pigmentosa and Stargardt disease. Expression of this gene is also associated with several types of cancer. This gene is expressed from at least five alternative promoters that are expressed in a tissue-dependent manner. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Mar 2009]

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of extracts from HT-29 cells using AMRe86843 at 1:1000.

注意事项 (Note)

For research use only .