

产品名称: Nanog Rabbit Monoclonal Antibody
产品货号: AMRe86847



产品概述 (Summary)

产品名称 (Production Name)	Nanog Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,FC,IP
种属反应性 (Reactivity)	Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Nanog
别名 (Alternative Names)	ENK; ecat4; 2410002E02Rik
基因 ID (Gene ID)	71950
蛋白 ID (SwissProt ID)	Q80Z64.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:200-1:500,ICC/IF 1:100-1:200,FC 1:50-1:100,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW:34 kDa; Observed MW:42 kDa

产品名称: Nanog Rabbit Monoclonal Antibody
产品货号: AMRe86847




研究背景 (Background)

The protein encoded by this gene is a DNA binding homeobox transcription factor involved in embryonic stem (ES) cell proliferation, renewal, and pluripotency. The encoded protein can block ES cell differentiation and can also autorepress its own expression in differentiating cells. Several transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Sep 2015]

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of extracts from F9 cells using AMRe86847 at 1:1000.

注意事项 (Note)

For research use only .