

产品名称: HES5 Rabbit Monoclonal Antibody
产品货号: AMRe86855



产品概述 (Summary)

产品名称 (Production Name)	HES5 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,IP
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	HES5
别名 (Alternative Names)	bHLHb38
基因 ID (Gene ID)	388585
蛋白 ID (SwissProt ID)	Q5TA89.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:5000,IHC 1:200-1:1000,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW:18 kDa; Observed MW:18 kDa

研究背景 (Background)


产品名称: HES5 Rabbit Monoclonal Antibody
产品货号: AMRe86855



This gene encodes a member of a family of basic helix-loop-helix transcriptional repressors. The protein product of this gene, which is activated downstream of the Notch pathway, regulates cell differentiation in multiple tissues. Disruptions in the normal expression of this gene have been associated with developmental diseases and cancer. [provided by RefSeq, Dec 2008]

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of extracts from Mouse small intestine tissue using AMRe86855 at 1:1000.

注意事项 (Note)

For research use only .