

产品名称: UBE2M Rabbit Monoclonal Antibody
产品货号: AMRe86857



产品概述 (Summary)

产品名称 (Production Name)	UBE2M Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,FC
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	UBE2M
别名 (Alternative Names)	UBC12; hUbc12; UBC-RS2
基因 ID (Gene ID)	9040
蛋白 ID (SwissProt ID)	P61081.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:2000-1:20000,IHC 1:100-1:200,ICC/IF 1:100-1:200,FC 1:100-1:1000
蛋白分子量 (Molecular Weight)	Calculated MW:21 kDa; Observed MW:21 kDa

研究背景 (Background)


产品名称: UBE2M Rabbit Monoclonal Antibody
产品货号: AMRe86857




The modification of proteins with ubiquitin is an important cellular mechanism for targeting abnormal or short-lived proteins for degradation. Ubiquitination involves at least three classes of enzymes: ubiquitin-activating enzymes, or E1s, ubiquitin-conjugating enzymes, or E2s, and ubiquitin-protein ligases, or E3s. This gene encodes a member of the E2 ubiquitin-conjugating enzyme family. The encoded protein is linked with a ubiquitin-like protein, NEDD8, which can be conjugated to cellular proteins, such as Cdc53/culin. [provided by RefSeq, Jul 2008]

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of extracts from Ramos cells using AMRe86857 at 1:1000.



Western blot analysis of extracts from HepG2, RAW264.7 cells and Rat kidney tissue using AMRe86857 at 1:5000.

注意事项 (Note)

For research use only .