

产品名称: CAPON Rabbit Monoclonal Antibody
产品货号: AMRe86861



产品概述 (Summary)

产品名称 (Production Name)	CAPON Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC/IF, FC
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	CAPON
别名 (Alternative Names)	CAPON; 6330408P19Rik
基因 ID (Gene ID)	9722
蛋白 ID (SwissProt ID)	O75052.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:5000, ICC/IF 1:20-1:50, FC 1:50-1:200
蛋白分子量 (Molecular Weight)	Calculated MW:56 kDa; Observed MW:56 kDa

研究背景 (Background)


产品名称: CAPON Rabbit Monoclonal Antibody
产品货号: AMRe86861



This gene encodes a cytosolic protein that binds to the signaling molecule, neuronal nitric oxide synthase (nNOS). This protein has a C-terminal PDZ-binding domain that mediates interactions with nNOS and an N-terminal phosphotyrosine binding (PTB) domain that binds to the small monomeric G protein, Dexras1. Studies of the related mouse and rat proteins have shown that this protein functions as an adapter protein linking nNOS to specific targets, such as Dexras1 and the synapsins. Alternative splicing results in multiple transcript variants encoding different isoforms. [provided by RefSeq, Sep 2009]

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of extracts from Rat brain cells using AMRe86861 at 1:1000.

注意事项 (Note)

For research use only .