

产品名称: UBE2D3 Rabbit Monoclonal Antibody
产品货号: AMRe86872



产品概述 (Summary)

产品名称 (Production Name)	UBE2D3 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IP
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	UBE2D3
别名 (Alternative Names)	UBC4/5; UBCH5C; E2(17)KB3
基因 ID (Gene ID)	7323
蛋白 ID (SwissProt ID)	P61077.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:5000,IP 1:10-1:100
蛋白分子量 (Molecular Weight)	Calculated MW:17 kDa; Observed MW:17 kDa

研究背景 (Background)


产品名称: UBE2D3 Rabbit Monoclonal Antibody
产品货号: AMRe86872




The modification of proteins with ubiquitin is an important cellular mechanism for targeting abnormal or short-lived proteins for degradation. Ubiquitination involves at least three classes of enzymes: ubiquitin-activating enzymes, or E1s, ubiquitin-conjugating enzymes, or E2s, and ubiquitin-protein ligases, or E3s. This gene encodes a member of the E2 ubiquitin-conjugating enzyme family. This enzyme functions in the ubiquitination of the tumor-suppressor protein p53, which is induced by an E3 ubiquitin-protein ligase. [provided by RefSeq, Jan 2017]

研究领域 (Research Area)

图片 (Image Data)



Western blot detection of UBE2D3 in K562, Rat Brain, C6, 3T3, HeLa cell lysates using UBE2D3 antibody (1:1000 diluted).



Western blot analysis of extracts from HepG2, RAW264.7 cells and Rat kidney tissue using AMRe86872 at 1:3000.

注意事项 (Note)

For research use only.