

产品名称: JAK2 Rabbit Monoclonal Antibody
产品货号: AMRe86890



产品概述 (Summary)

产品名称 (Production Name)	JAK2 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC/IF, IP
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	JAK2
别名 (Alternative Names)	JTK10; THCYT3
基因 ID (Gene ID)	3717
蛋白 ID (SwissProt ID)	O60674.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:5000, ICC/IF 1:100-1:200, IP 1:50-1:100
蛋白分子量 (Molecular Weight)	Calculated MW:131 kDa; Observed MW:131 kDa

研究背景 (Background)


产品名称: JAK2 Rabbit Monoclonal Antibody
产品货号: AMRe86890




This gene product is a protein tyrosine kinase involved in a specific subset of cytokine receptor signaling pathways. It has been found to be constitutively associated with the prolactin receptor and is required for responses to gamma interferon. Mice that do not express an active protein for this gene exhibit embryonic lethality associated with the absence of definitive erythropoiesis. [provided by RefSeq, Jul 2008]

研究领域 (Research Area)

图片 (Image Data)



Western blot detection of JAK2 in K562,3T3 cell lysates using JAK2 antibody(1:1000 diluted).



Western blot analysis of extracts from Rat spleen tissue using AMRe86890 at 1:1000.

注意事项 (Note)

For research use only .