

产品名称: Proteasome 26S S10 Rabbit Monoclonal Antibody
产品货号: AMRe86893



产品概述 (Summary)

产品名称 (Production Name)	Proteasome 26S S10 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,FC
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Proteasome 26S S10
别名 (Alternative Names)	p28; p28(GANK); dJ889N15.2
基因 ID (Gene ID)	5716
蛋白 ID (SwissProt ID)	O75832.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:5000,IHC 1:50-1:100,ICC/IF 1:100-1:200,FC 1:50-1:100
蛋白分子量 (Molecular Weight)	Calculated MW:24 kDa; Observed MW:24 kDa

研究背景 (Background)


产品名称: Proteasome 26S S10 Rabbit Monoclonal Antibody
产品货号: AMRe86893



This gene encodes a subunit of the PA700/19S complex, which is the regulatory component of the 26S proteasome. The 26S proteasome complex is required for ubiquitin-dependent protein degradation. This protein is a non-ATPase subunit that may be involved in protein-protein interactions. Aberrant expression of this gene may play a role in tumorigenesis. Two transcripts encoding different isoforms have been described. Pseudogenes have been identified on chromosomes 3 and 20.[provided by RefSeq, Mar 2011]

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of extracts from HeLa,SW620 cells using AMRe86893 at 1:4000.

注意事项 (Note)

For research use only .