

产品名称: TAOK1 Rabbit Monoclonal Antibody
产品货号: AMRe86894



产品概述 (Summary)

产品名称 (Production Name)	TAOK1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	TAOK1
别名 (Alternative Names)	PSK2; TAO1; KFC-B; MARKK; PSK-2; hKFC-B; hTAOK1; MAP3K16
基因 ID (Gene ID)	57551
蛋白 ID (SwissProt ID)	Q7L7X3.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000
蛋白分子量 (Molecular Weight)	Calculated MW:116 kDa; Observed MW:116 kDa

研究背景 (Background)


产品名称: TAOK1 Rabbit Monoclonal Antibody
产品货号: AMRe86894




Serine/threonine-protein kinase involved in various processes such as p38/MAPK14 stress-activated MAPK cascade, DNA damage response and regulation of cytoskeleton stability

研究领域 (Research Area)

图片 (Image Data)



Western blot detection of TAOK1 in C6,Ramos using TAOK1 antibody(1:1000 diluted)



Western blot analysis of extracts from HepG2, RAW264.7 cells and Rat spleen tissue using AMRe86894 at 1:1000.

注意事项 (Note)

For research use only .