

产品名称: Tin2 Rabbit Monoclonal Antibody
产品货号: AMRe86896



产品概述 (Summary)

产品名称 (Production Name)	Tin2 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Tin2
别名 (Alternative Names)	TIN2; DKCA3
基因 ID (Gene ID)	26277
蛋白 ID (SwissProt ID)	Q9BSI4.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000
蛋白分子量 (Molecular Weight)	Calculated MW:50 kDa; Observed MW:39 kDa

研究背景 (Background)


产品名称: Tin2 Rabbit Monoclonal Antibody
产品货号: AMRe86896




This gene encodes one of the proteins of the shelterin, or telosome, complex which protects telomeres by allowing the cell to distinguish between telomeres and regions of DNA damage. The protein encoded by this gene is a critical part of shelterin; it interacts with the three DNA-binding proteins of the shelterin complex, and it is important for assembly of the complex. Mutations in this gene cause dyskeratosis congenita (DKC), an inherited bone marrow failure syndrome. [provided by RefSeq, Mar 2010]

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of extracts from MCF-7, HUVEC cells using AMRe86896 at 1:1000.



Western blot analysis of extracts from C2C12, T24 cells and Rat spleen tissue using AMRe86896 at 1:1000.

注意事项 (Note)

For research use only .