

产品名称: CDCA5 Rabbit Monoclonal Antibody
产品货号: AMRe86898



产品概述 (Summary)

产品名称 (Production Name)	CDCA5 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,FC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	CDCA5
别名 (Alternative Names)	SORORIN
基因 ID (Gene ID)	113130
蛋白 ID (SwissProt ID)	Q96FF9.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:5000,IHC 1:50-1:100,ICC/IF 1:100-1:200,FC 1:50-1:200
蛋白分子量 (Molecular Weight)	Calculated MW:28 kDa; Observed MW:35 kDa

研究背景 (Background)


产品名称: CDCA5 Rabbit Monoclonal Antibody
产品货号: AMRe86898




Regulator of sister chromatid cohesion in mitosis stabilizing cohesin complex association with chromatin. May antagonize the action of WAPL which stimulates cohesin dissociation from chromatin. Cohesion ensures that chromosome partitioning is accurate in both meiotic and mitotic cells and plays an important role in DNA repair. Required for efficient DNA double-stranded break repair.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of extracts from HeLa cells using AMRe86898 at 1:1000.



Western blot analysis of extracts from C2C12, T24 cells using AMRe86898 at 1:3000.

注意事项 (Note)

For research use only .