

产品名称: SRA1 Rabbit Monoclonal Antibody
产品货号: AMRe86904



产品概述 (Summary)

产品名称 (Production Name)	SRA1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,FC
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	SRA1
别名 (Alternative Names)	SRA; SRAP; STRAA1; pp7684
基因 ID (Gene ID)	10011
蛋白 ID (SwissProt ID)	Q9HD15.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:5000,IHC 1:50-1:100,ICC/IF 1:200-1:500,FC 1:100-1:200
蛋白分子量 (Molecular Weight)	Calculated MW:26 kDa; Observed MW:36 kDa

研究背景 (Background)


产品名称: SRA1 Rabbit Monoclonal Antibody
产品货号: AMRe86904



Both long non-coding and protein-coding RNAs are transcribed from this gene, and they represent alternatively spliced transcript variants. This gene was initially defined as a non-coding RNA, which is a coactivator for several nuclear receptors (NRs) and is associated with breast cancer. It has now been found that this gene is involved in the regulation of many NR and non-NR activities, including metabolism, adipogenesis and chromatin organization. The long non-coding RNA transcripts interact with a variety of proteins, including the protein encoded by this gene. The encoded protein acts as a transcriptional repressor by binding to the non-coding RNA. [provided by RefSeq, Mar 2012]

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of extracts from Rat liver tissue using AMRe86904 at 1:1000.

注意事项 (Note)

For research use only .