

产品名称: GBA Rabbit Monoclonal Antibody
产品货号: AMRe86910



产品概述 (Summary)

产品名称 (Production Name)	GBA Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	GBA
别名 (Alternative Names)	GCB; GBA1; GLUC
基因 ID (Gene ID)	2629
蛋白 ID (SwissProt ID)	P04062.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:5000,IHC 1:50-1:100
蛋白分子量 (Molecular Weight)	Calculated MW:60 kDa; Observed MW:60 kDa

研究背景 (Background)


产品名称: GBA Rabbit Monoclonal Antibody
产品货号: AMRe86910



This gene encodes a lysosomal membrane protein that cleaves the beta-glucosidic linkage of glycosylceramide, an intermediate in glycolipid metabolism. Mutations in this gene cause Gaucher disease, a lysosomal storage disease characterized by an accumulation of glucocerebrosides. A related pseudogene is approximately 12 kb downstream of this gene on chromosome 1. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jan 2010]

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of extracts from MCF-7 cells using AMRe86910 at 1:1000.

注意事项 (Note)

For research use only .