

产品名称: LMCD1 Rabbit Monoclonal Antibody
产品货号: AMRe86913



产品概述 (Summary)

产品名称 (Production Name)	LMCD1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC/IF, IP
种属反应性 (Reactivity)	Human, Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	LMCD1
别名 (Alternative Names)	Dyxin; LIM and cysteine-rich domains protein 1
基因 ID (Gene ID)	29995
蛋白 ID (SwissProt ID)	Q9NZU5.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:5000, ICC/IF 1:100-1:200, IP 1:10-1:100
蛋白分子量 (Molecular Weight)	Calculated MW:41 kDa; Observed MW:

研究背景 (Background)


产品名称: LMCD1 Rabbit Monoclonal Antibody
产品货号: AMRe86913



This gene encodes a member of the LIM-domain family of zinc finger proteins. The encoded protein contains an N-terminal cysteine-rich domain and two C-terminal LIM domains. The presence of LIM domains suggests involvement in protein-protein interactions. The protein may act as a co-regulator of transcription along with other transcription factors. Alternate splicing results in multiple transcript variants of this gene. [provided by RefSeq, May 2013]

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of extracts from C2C12 cells using AMRe86913 at 1:1000.

注意事项 (Note)

For research use only .