

产品名称: PPAR gamma Rabbit Monoclonal Antibody
产品货号: AMRe86915



产品概述 (Summary)

产品名称 (Production Name)	PPAR gamma Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC/IF, FC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	PPAR gamma
别名 (Alternative Names)	GLM1; CIMT1; NR1C3; PPARG1; PPARG2; PPARGgamma
基因 ID (Gene ID)	5468
蛋白 ID (SwissProt ID)	P37231.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000, ICC/IF 1:100-1:200, FC 1:200-1:1000
蛋白分子量 (Molecular Weight)	Calculated MW:58 kDa; Observed MW:58 kDa

研究背景 (Background)


产品名称: PPAR gamma Rabbit Monoclonal Antibody
产品货号: AMRe86915




This gene encodes a member of the peroxisome proliferator-activated receptor (PPAR) subfamily of nuclear receptors. PPARs form heterodimers with retinoid X receptors (RXRs) and these heterodimers regulate transcription of various genes. Three subtypes of PPARs are known: PPAR-alpha, PPAR-delta, and PPAR-gamma. The protein encoded by this gene is PPAR-gamma and is a regulator of adipocyte differentiation. Additionally, PPAR-gamma has been implicated in the pathology of numerous diseases including obesity, diabetes, atherosclerosis and cancer. Alternatively spliced transcript variants that encode different isoforms have been described. [provided by RefSeq, Jul 2008]

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of extracts from Human fetal tissue using AMRe86915 at 1:1000.



Western blot analysis of extracts from K562 cells using AMRe86915 at 1:1000.

注意事项 (Note)

For research use only .