

产品名称: Lysophospholipase 1 Rabbit Monoclonal Antibody
产品货号: AMRe86927



产品概述 (Summary)

产品名称 (Production Name)	Lysophospholipase 1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Lysophospholipase 1
别名 (Alternative Names)	APT1; LPL1; APT-1; LPL-I; hAPT1
基因 ID (Gene ID)	10434
蛋白 ID (SwissProt ID)	O75608.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:2000-1:20000
蛋白分子量 (Molecular Weight)	Calculated MW:25 kDa; Observed MW:25 kDa

研究背景 (Background)


产品名称: Lysophospholipase 1 Rabbit Monoclonal Antibody
产品货号: AMRe86927




This gene encodes a member of the alpha/beta hydrolase superfamily. The encoded protein functions as a homodimer, exhibiting both depalmitoylating as well as lysophospholipase activity, and may be involved in Ras localization and signaling. Alternate splicing results in multiple transcript variants. Pseudogenes of this gene have been defined on chromosomes 4, 6, and 7. [provided by RefSeq, Jul 2013]

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of extracts from Mouse spleen tissue using AMRe86927 at 1:1000.



Western blot analysis of extracts from RAW264.7 cells using AMRe86927 at 1:5000.

注意事项 (Note)

For research use only .