

产品名称: LAMP2A Rabbit Monoclonal Antibody
产品货号: AMRe86977



产品概述 (Summary)

产品名称 (Production Name)	LAMP2A Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,FC,IP
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	LAMP2A
别名 (Alternative Names)	LAMPB; CD107b; LAMP-2; LGP-96; LGP110
基因 ID (Gene ID)	3920
蛋白 ID (SwissProt ID)	P13473.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:20-1:100,FC 1:100-1:200,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW:45 kDa; Observed MW:110 kDa

研究背景 (Background)


产品名称: LAMP2A Rabbit Monoclonal Antibody
产品货号: AMRe86977




The protein encoded by this gene is a member of a family of membrane glycoproteins. This glycoprotein provides selectins with carbohydrate ligands. It may play a role in tumor cell metastasis. It may also function in the protection, maintenance, and adhesion of the lysosome. Alternative splicing of this gene results in multiple transcript variants encoding distinct proteins. [provided by RefSeq, Jul 2008]

研究领域 (Research Area)


图片 (Image Data)



Western blot detection of LAMP2A in 3T3 cell lysates using LAMP2A antibody(1:1000 diluted).



Western blot analysis of extracts from HeLa, MCF-7 cells using AMRe86977 at 1:1000.



产品名称: LAMP2A Rabbit Monoclonal Antibody
产品货号: AMRe86977



Immunohistochemical analysis of paraffin-embedded human lung cancer using AMRe86977 antibody.

注意事项 (Note)

For research use only .