

产品名称: IKK alpha Rabbit Monoclonal Antibody
产品货号: AMRe86988



产品概述 (Summary)

产品名称 (Production Name)	IKK alpha Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,FC,IP
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	IKK alpha
别名 (Alternative Names)	IKK1; IKKA; IKBKA; TCF16; NFKBIKA; IKK-alpha
基因 ID (Gene ID)	1147
蛋白 ID (SwissProt ID)	O15111.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:50-1:100,FC 1:20-1:50,IP 1:50-1:100
蛋白分子量 (Molecular Weight)	Calculated MW:85 kDa; Observed MW:85 kDa

研究背景 (Background)


产品名称: IKK alpha Rabbit Monoclonal Antibody
产品货号: AMRe86988




This gene encodes a member of the serine/threonine protein kinase family. The encoded protein, a component of a cytokine-activated protein complex that is an inhibitor of the essential transcription factor NF-kappa-B complex, phosphorylates sites that trigger the degradation of the inhibitor via the ubiquination pathway, thereby activating the transcription factor. [provided by RefSeq, Jul 2008]

研究领域 (Research Area)

图片 (Image Data)



Western blot detection of IKK alpha in K562,C6,3T3,HeLa cell lysates using IKK alpha antibody(1:1000 diluted).



Western blot analysis of extracts from Mouse kidney tissue using AMRe86988 at 1:1000.

注意事项 (Note)

For research use only .