

产品名称: FNDC5 Rabbit Monoclonal Antibody
产品货号: AMRe86999



产品概述 (Summary)

产品名称 (Production Name)	FNDC5 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	FNDC5
别名 (Alternative Names)	FRCP2; irisin
基因 ID (Gene ID)	252995
蛋白 ID (SwissProt ID)	Q8NAU1.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:5000
蛋白分子量 (Molecular Weight)	Calculated MW:24 kDa; Observed MW:21 kDa

研究背景 (Background)


产品名称: FNDC5 Rabbit Monoclonal Antibody
产品货号: AMRe86999




This gene encodes a secreted protein that is released from muscle cells during exercise. The encoded protein may participate in the development of brown fat. Translation of the precursor protein initiates at a non-AUG start codon at a position that is conserved as an AUG start codon in other organisms. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jun 2013]

研究领域 (Research Area)

图片 (Image Data)



Western blot detection of FNDC5 in K562, C6, 3T3, HeLa cell lysates using FNDC5 antibody (1:1000 diluted).



Western blot analysis of extracts from Mouse muscle tissue using AMRe86999 at 1:1000.

注意事项 (Note)

For research use only .