

产品名称: RIP Rabbit Monoclonal Antibody
产品货号: AMRe87001



产品概述 (Summary)

产品名称 (Production Name)	RIP Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,FC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	RIP
别名 (Alternative Names)	RIP; RIP1; RIP-1
基因 ID (Gene ID)	8737
蛋白 ID (SwissProt ID)	Q13546.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:5000,FC 1:10-1:100
蛋白分子量 (Molecular Weight)	Calculated MW:76 kDa; Observed MW:76 kDa

研究背景 (Background)


产品名称: RIP Rabbit Monoclonal Antibody
产品货号: AMRe87001



This gene encodes a member of the receptor-interacting protein (RIP) family of serine/threonine protein kinases. The encoded protein plays a role in inflammation and cell death in response to tissue damage, pathogen recognition, and as part of developmental regulation. RIPK1/RIPK3 kinase-mediated necrosis is referred to as necroptosis. Genetic disruption of this gene in mice results in death shortly after birth. [provided by RefSeq, Aug 2017]

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of extracts from HepG2, T24 cells using AMRe87001 at 1:2000.

注意事项 (Note)

For research use only .