

产品名称: GRB7 Rabbit Monoclonal Antibody
产品货号: AMRe87006



产品概述 (Summary)

产品名称 (Production Name)	GRB7 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	GRB7
别名 (Alternative Names)	B47; Epidermal growth factor receptor GRB-7; GRB7 adapter protein
基因 ID (Gene ID)	2886
蛋白 ID (SwissProt ID)	Q14451.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:10000,IHC 1:200-1:500
蛋白分子量 (Molecular Weight)	Calculated MW:60 kDa; Observed MW:60 kDa

研究背景 (Background)


产品名称: GRB7 Rabbit Monoclonal Antibody
产品货号: AMRe87006



The product of this gene belongs to a small family of adapter proteins that are known to interact with a number of receptor tyrosine kinases and signaling molecules. This gene encodes a growth factor receptor-binding protein that interacts with epidermal growth factor receptor (EGFR) and ephrin receptors. The protein plays a role in the integrin signaling pathway and cell migration by binding with focal adhesion kinase (FAK). Several transcript variants encoding two different isoforms have been found for this gene. [provided by RefSeq, Jun 2011]

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of extracts from HepG2 cells using AMRe87006 at 1:1000.

注意事项 (Note)

For research use only .