

产品名称: BAP1 Rabbit Monoclonal Antibody
产品货号: AMRe87019



产品概述 (Summary)

产品名称 (Production Name)	BAP1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,FC
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	BAP1
别名 (Alternative Names)	UCLH2; hucep-6; HUCEP-13
基因 ID (Gene ID)	8314
蛋白 ID (SwissProt ID)	Q92560.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:50-1:200,FC 1:50-1:100
蛋白分子量 (Molecular Weight)	Calculated MW:80 kDa; Observed MW:95 kDa

研究背景 (Background)


产品名称: BAP1 Rabbit Monoclonal Antibody
产品货号: AMRe87019




This gene belongs to the ubiquitin C-terminal hydrolase subfamily of deubiquitinating enzymes that are involved in the removal of ubiquitin from proteins. The encoded enzyme binds to the breast cancer type 1 susceptibility protein (BRCA1) via the RING finger domain of the latter and acts as a tumor suppressor. In addition, the enzyme may be involved in regulation of transcription, regulation of cell cycle and growth, response to DNA damage and chromatin dynamics. Germline mutations in this gene may be associated with tumor predisposition syndrome (TPDS), which involves increased risk of cancers including malignant mesothelioma, uveal melanoma and cutaneous melanoma. [provided by RefSeq, May 2013]

研究领域 (Research Area)

图片 (Image Data)



Western blot detection of BAP1 in HeLa, A549, HL-60, U2OS cell lysates using BAP1 antibody (1:1000 diluted).



Western blot analysis of extracts from A431 cells using AMRe87019 at 1:1000.

注意事项 (Note)

For research use only .